



# S R ENDOSCOPY SYSTEM

Otolaryngology Product







## S R ENDOSCOPY SYSTEM

Otolaryngology Product

WE INTRODUCE OUR SELF AS MEDICAL EQUIPMENT COMPANY, WE DEALS MOSTLY IN ENT EQUIPMENTS, LAPAROSCOPY EQUIPMENTS ORTHO, AND DENTAL MICROSCOPES. ALL THESE PRODUCTS ARE MADE BY OUR SELVES. SO THAT WE CAN PROVIDE OUR CUSTOMERS WITH GOOD EQUIPMENT AND GOOD SERVICE, OUR OFFICE IS LOCATED IN HARIDWAR. WHICH COMES UNDER THE STATE OF UTTARAKHAND. WE ALSO SPECIALIZE IN REPAIRING SINOSCOPE, LAPROSCOPE, CYSTOSCOPE OF ANY COMPANY WE ALWAYS RESPECT OUR CUSTOMERS & WILL YOUR TIME IS VALUABLE FOR US. PLEASE GIVE A CHANCE OF SERVICE. WELCOME TO CALL US FOR INFORMATION OR ANY KIND OF MEDICAL EQUIPMENT.

 Mantra Society Flat No. D-708, Near Raja Biscuit Chowk, Opp - Sec-2  
Sidcul Haridwar -249403 (U.K)

 09639567270, 08279491608

Email : [srendoscopysystem@gmail.com](mailto:srendoscopysystem@gmail.com)



## **COMPACT ENDOSCOPY UNIT INCLUDING OF ENDOSCOPY CAMERA S-1 L.E.D PORTABLE LIGHT SOURCE 8 WATT L.E.D MONITOR 11 INCH**

- Easy to Carry for Other Place
- Start with single Power Connection
- Double video out Put for extra monitor
- Video System (CVBS))
- Images Sencor 1/3" CCD
- Auto White Balance
- Weight Approx 4 Kg



## **Endoscopy Camera Model No. S-1**

- Useful For Diagnostic Procedure
- Single Video Out Put
- With Economical Price



## **Endoscopy Camera (Standrad) Model No. S-2 For Diagnostic & Surgical Procedure**

- Video System (CVBS)
- Resolution 700 TVL
- Double Video Out for Extra Monitor
- Auto White Balance
- Power 100-240 Vac





## Endoscopy Camera Mode No. S-1

- Useful For Diagnostic Procedure
- Single Video OutPut
- With economical Price



## Refurbished Cameras Single & Three Chip HD Storz & Stryker



## LED Light Source 100/150 watt

- Luminous Flux -4500 To 6000
- Bright & White Light
- L.E.D Life 50000 Hours
- Power Supply 110V-220V
- Fiber optic cable Fitting Storz Type



- Halogen Light Source 250 watt With Twin Bulb






## Portable USB L.E.D Light Source 10 Watt

- Can be connect any kind of Sinoscope
- Non Rechargeable
- Direct Connect to Power Point
- Can be use with any Power Bank & Laptop
- L.E.D Life approx 30000 Hours.



## ENT MICRO DIBRIDER

- Xomed Fitting Blade
- Detachable Hand Piece
- Clock Wise, Anticlockwise & Oscillating Mode
- Very Light Weighted Hand Piece
- Smooth Cutting



## Cautery Machine 250/400 Watt

- Monopolar Function and Bipolar Function Available
- Pure and blend Cut function
- Useful For all kind of Surgery



## Fiber optic Cable

- Assembled with top Quality Fibers
- Cable Length 7.5 feet  
Diameter from 3.5 to 5mm Available





## Portable L.E.D Endoscopy Light Source

- Can be Connect any kind of Sinoscope
- Rechargeable Battery
- L.E.D Life Approx 30000 Hours.

## ENT Operative Microscope Specification

- Binocular 90<sup>0</sup> Stright
- Magnification Charges 3 Step
- Total Magnification 5x10x20x
- Working Distance 200mm
- Objective Lens f= 200
- Heavy Floor Stand With caster Wheel
- With L.E.D. Light Source



## Sinoscope, Otoenloscope, Laryngoscope

- 0 Degree 4mm
- 0 Degree 2.7mm
- 0 Degree Otoendoscope
- 30 Degree 4mm Sinoscope
- 45 Degree 4mm Sinoscope
- 70 Degree 4mm Sinoscope
- 90 Degree 8mm Laryngoscope



## L.E.D Soft Head Light rechargeable (8 Watt)





## ENT OPD Unit

- Inbuilt light Source 15v/ 150 watt
- Inbuilt Suction Machine with twin jars
- Xray View Box
- Nasal Cautery with Temperature Control
- Mirror Warmer
- Instrument Trey with cover 2 Nos
- Optical Fiber Cable
- Head band for fiber cable Attachment
- Sliding Drawer for Storage
- Hanger for head Light
- Sinoscope Holder
- Otoscope Holder
- Container for waste gojj Piece



## Endoscopy Camera Model No. S-3

- Useful for Operative Procedure
- White Balance in Camera Console
- Resolution 1000 TVL
- Double Video Out Put



- ENT Micro Motor With NSK Hand Piece





**S R ENDOSCOPY SYSTEM**

# **S R ENDOSCOPY SYSTEM**




**Repairing : IMP : We repair all make of Sinoscope, Laproscope, Cystoscope etc.**

Facility for repair of all kind of Endoscopy equipment Which includes, Camera, Telescope, Microscope, Debrider, Flexible Gastroscope , CO2 Insuflator, Light Source, Fiber Optic Cable, Hand instruments etc.



 Mantra Society Flat No. D-708,  
Near Raja Biscuit Chowk,  
Opp - Sec-2, Sidcul  
Haridwar -249403 (U.K)

 09639567270, 08279491608

 [srendoscopysystem@gmail.com](mailto:srendoscopysystem@gmail.com)